


Certificate of Analysis

(C of A)

Name:		ATFB, SE 4-Azido-2,3,5,6-tetrafluorobenzoic acid, succinimidyl ester						
Catalog Number:		6977	Lot:		50829			
CAS Number:		126695-58-7	Formula:		C ₁₁ H ₈ F ₄ N ₄ O ₄			
Molecular Weight:		332.17 (g/mol)						
		Observed		Standard				
HPLC Purity:		97% at 254 nm		≥ 90% at 254 nm				
Absorption (Ethanol)								
Wavelength (λ _{max}):		274 nm		274 ± 3 nm				
Extinction:		25,900 cm ⁻¹ M ⁻¹		24,000 ± 2,500 cm ⁻¹ M ⁻¹				
TLC Purity:		> 99%		≥ 90%				
TLC System:		2:1 Hexane / Ethyl acetate						
Storage:		Freeze: Store at -5 to -30°C (F), Protect from moisture (D), Protect from light (L)						
Production Lot:		14354	Release:		27-Sep-2019	Expiration:		15-Jul-2019
Structure:								

This data reflects the accuracy of the analysis of the subject lot
 For research use only. Not for human use. No animal source.
 This material is chemically synthesized and manufactured in the USA

Approved by: Dianna M Wineinger
 Analytical Services Manager

Date: 27-Sep-2019

Signature: Dianna M Wineinger