


Certificate of Analysis (C of A)					
Name:	D-Biotin				
Catalog Number:	6738		Lot:	51144	
CAS Number:	58-85-5		Formula:	C ₁₀ H ₁₆ N ₂ O ₃ S	
Molecular Weight:	244.31 (g/mol)				
	Observed		Standard		
TLC Purity:	> 99%		≥ 98%		
TLC System:	80:20 Chloroform / Methanol (with Biotin dip)				
Storage:	Store at Room Temperature (RT), Protect from moisture (D), Protect from light (L)				
Production Lot:	15761	Release:	20-Apr-2017	Expiration:	20-Apr-2037
Structure:					

This data reflects the accuracy of the analysis of the subject lot
 For research use only. Not for human use. No animal source.
 This material is chemically synthesized and manufactured in the USA

Approved by: Dianna M Wineinger
 Analytical Services Manager

Date: 20-Apr-2017

Signature: 