


Certificate of Analysis (C of A)					
Name:	N-(2-Aminoethyl)maleimide, trifluoroacetic acid salt C2-Maleimide, TFA				
Catalog Number:	7018	Lot:	51395		
CAS Number:	146474-00-2	Formula:	C ₈ H ₉ F ₃ N ₂ O ₄		
Molecular Weight:	254.16 (g/mol)				
	Observed	Standard			
TLC Purity:	> 98%		≥ 90%		
TLC System:	60:30:7:3 Chloroform / Methanol / Water / Acetic Acid with Ninhydrin Dip				
Storage:	Freeze: Store at -5 to -30°C (F), Protect from moisture (D), Protect from light (L)				
Production Lot:	15265	Release:	15-Jan-2016	Expiration:	15-Jan-2026
Structure:					

This data reflects the accuracy of the analysis of the subject lot
 For research use only. Not for human use. No animal source.
 This material is chemically synthesized and manufactured in the USA

Approved by: Dianna M Wineinger
 Analytical Services Manager

Date: 15-Jan-2016

Signature: 