


Certificate of Analysis					
(C of A)					
Name:		SBFI, AM *Sp-Pak* *Cell Permeant*			
Catalog Number:		6212	Lot:		52325
CAS Number:		129423-53-6	Formula:		C ₅₆ H ₅₈ N ₂ O ₂₃
Molecular Weight:		1127.03 (g/mol)			
	Observed		Standard		
HPLC Purity:		99% at 254 nm	≥ 95% at 254 nm		
Absorption (Methanol)					
Wavelength (λ _{max}):		379 nm	380 ± 3 nm		
Extinction:		Set 34,000 cm ⁻¹ M ⁻¹	Set 34,000 ± 3,500 cm ⁻¹ M ⁻¹		
TLC Purity:		> 99%	≥ 95%		
TLC System:		90:10 Chloroform / Methanol			
Storage:		Freeze: Store at -5 to -30°C (F), Protect from moisture (D), Protect from light (L)			
Production Lot:		15533	Release:		19-Jan-2023
			Expiration:		19-Jan-2033
Structure:					

This data reflects the accuracy of the analysis of the subject lot
 For research use only. Not for human use. No animal source.
 This material is chemically synthesized and manufactured in the USA

Approved by: Nabi
 Analytical Services Manager
Signature: 

Date: 19-Jan-2023