


Certificate of Analysis (C of A)					
Name:	CNF, SE *Mixed Isomers* 5(6)-Carboxynaphthofluorescein , succinimidyl ester				
Catalog Number:	6964		Lot:	52356	
CAS Number:	150347-58-3		Formula:	C ₃₃ H ₁₉ NO ₉	
Molecular Weight:	573.51 (g/mol)				
	Observed		Standard		
HPLC Purity:	92% at 254 nm		≥ 90% at 254 nm		
Absorption (Buffer pH10)					
Wavelength (λ _{max}):	602 nm		602 ± 3 nm		
Extinction:	45,800 cm ⁻¹ M ⁻¹		48,000 ± 4,000 cm ⁻¹ M ⁻¹		
TLC Purity:	> 98%		≥ 90%		
TLC System:	90:10 Chloroform / Methanol				
Storage:	Freeze: Store at -5 to -30°C (F), Protect from moisture (D), Protect from light (L)				
Production Lot:	11937	Release:	01-Jul-2014	Expiration:	03-Apr-2028
Structure:					

This data reflects the accuracy of the analysis of the subject lot
 For research use only. Not for human use. No animal source.
 This material is chemically synthesized and manufactured in the USA

Approved by: Dianna M Wineinger
 Analytical Services Manager

Date: 01-Jul-2014

Signature: 