



Certificate of Analysis

(C of A)

Name:		MESG 7-Methyl-6-thioguanosine				
Catalog Number:		6764	Lot:		52387	
CAS Number:		55727-10-1	Formula:		C ₁₁ H ₁₅ N ₅ O ₄ S	
Molecular Weight:		313.33 (g/mol)				
		Observed		Standard		
HPLC Purity:		93% at 254 nm		≥ 85% at 254 nm		
Absorption (Buffer pH 7)						
Wavelength (λ _{max}):		332 nm		333 ± 3 nm		
Extinction:		16,600 cm ⁻¹ M ⁻¹		17,000 ± 1,500 cm ⁻¹ M ⁻¹		
TLC Purity:		> 95%		≥ 90%		
TLC System:		80:20 Acetonitrile / Water				
Storage:		Freeze: Store at -5 to -30°C (F), Protect from moisture (D), Protect from light (L)				
Production Lot:		16245	Release:		14-Dec-2023	
				Expiration:		14-Dec-2033
Structure:						

This data reflects the accuracy of the analysis of the subject lot
 For research use only. Not for human use. No animal source.
 This material is chemically synthesized and manufactured in the USA

Approved by: Nabi
 Analytical Services Manager
Signature: 

Date: 14-Dec-2023