


Certificate of Analysis (C of A)					
Name:	AMBH 2-Acetamido-4-mercaptobutanoic acid hydrazide				
Catalog Number:	6470	Lot:	51864		
CAS Number:	77076-41-6	Formula:	C ₆ H ₁₃ N ₃ O ₂ S		
Molecular Weight:	191.25 (g/mol)				
	Observed		Standard		
TLC Purity:	> 99%		≥ 98%		
TLC System:	75:20:5 Chloroform/ Methanol/ Acetic acid				
Storage:	Freeze: Store at -5 to -30°C (F), Protect from moisture (D), Protect from light (L)				
Production Lot:	11667	Release:	29-Oct-2014	Expiration:	23-Mar-2023
Structure:	 <chem>CC(=O)NCCSCCNC(=O)N</chem>				

This data reflects the accuracy of the analysis of the subject lot
 For research use only. Not for human use. No animal source.
 This material is chemically synthesized and manufactured in the USA

Approved by: Dianna M Wineinger
 Analytical Services Manager

Date: 29-Oct-2014

Signature: 