



| Certificate of Analysis (C of A) | | | | | |
|-------------------------------------|---|----------|-------------|---|-------------|
| Name: | CNF, SE *Mixed Isomers* 5(6)-Carboxynaphthofluorescein , succinimidyl ester | | | | |
| Catalog Number: | 6964 | | Lot: | 51708 | |
| CAS Number: | 150347-58-3 | | Formula: | C ₃₃ H ₁₉ NO ₉ | |
| Molecular Weight: | 573.51 (g/mol) | | | | |
| | | | | | |
| | Observed | | | Standard | |
| HPLC Purity: | 94% at 254 nm | | | ≥ 90% at 254 nm | |
| Absorption (Buffer pH10) | | | | | |
| Wavelength (λ _{max}): | 602 nm | | | 602 ± 3 nm | |
| Extinction: | 47,500 cm ⁻¹ M ⁻¹ | | | 48,000 ± 4,000 cm ⁻¹ M ⁻¹ | |
| TLC Purity: | > 98% | | | ≥ 90% | |
| TLC System: | 90:10 Chloroform / Methanol | | | | |
| Storage: | Freeze: Store at -5 to -30°C (F), Protect from moisture (D), Protect from light (L) | | | | |
| Production Lot: | 16256 | Release: | 23-Feb-2018 | Expiration: | 23-Feb-2023 |
| Structure: |  | | | | |

This data reflects the accuracy of the analysis of the subject lot
 For research use only. Not for human use. No animal source.
 This material is chemically synthesized and manufactured in the USA

Approved by: Dianna M Wineinger
 Analytical Services Manager

Date: 23-Feb-2018

Signature: 

Date: 09-Dec-2022