


Certificate of Analysis (C of A)			
Name:	FDGlu Fluorescein di-β-D-glucopyranoside		
Catalog Number:	6912	Lot:	51533
CAS Number:	129787-66-2	Formula:	C ₃₂ H ₃₂ O ₁₅
Molecular Weight:	656.60 (g/mol)		
	Observed	Standard	
HPLC Purity:	96% at 254 nm		≥ 90% at 254 nm
Absorption (Methanol)			
Wavelength (λ_{max}):	273 nm		273 ± 3 nm
Extinction:	6,400 cm ⁻¹ M ⁻¹		6,500 ± 600 cm ⁻¹ M ⁻¹
Storage:	Freeze: Store at -5 to -30°C (F), Protect from moisture (D), Protect from light (L)		
Production Lot:	15840	Release:	10-Jul-2017 Expiration: 10-Jul-2027
Structure:			

This data reflects the accuracy of the analysis of the subject lot
 For research use only. Not for human use. No animal source.
 This material is chemically synthesized and manufactured in the USA

Approved by: Dianna M Wineinger
 Analytical Services Manager

Date: 10-Jul-2017

Signature:

