


Certificate of Analysis (C of A)					
Name:	Octadecyl 4-chlorobenzenesulfonate Octadecyl p-chlorobenzenesulfonate				
Catalog Number:	7058		Lot:	50677	
CAS Number:	34184-41-3		Formula:	C ₂₄ H ₄₁ ClO ₃ S	
Molecular Weight:	445.10 (g/mol)				
	Observed		Standard		
HPLC Purity:	99% at 254 nm		≥ 95% at 254 nm		
TLC Purity:	> 99%		≥ 95%		
TLC System:	1:1 Hexane/ Chloroform + Phosphomolybdic Acid (PMA)				
Melting Point:	63-65 C		64-66 C		
Storage:	Store at Room Temperature (RT), Protect from moisture (D), Protect from light (L)				
Production Lot:	14720	Release:	20-Jan-2015	Expiration:	20-Jan-2035
Structure:					

This data reflects the accuracy of the analysis of the subject lot
 For research use only. Not for human use. No animal source.
 This material is chemically synthesized and manufactured in the USA

Approved by: Dianna M Wineinger
 Analytical Services Manager

Date: 20-Jan-2015

Signature: 