


Certificate of Analysis (C of A)					
Name:	CellHunt Blue CMHC 7-Hydroxy-4-chloromethylcoumarin				
Catalog Number:	7135		Lot:	52165	
CAS Number:	25392-41-0		Formula:	C ₁₀ H ₇ ClO ₃	
Molecular Weight:	210.61 (g/mol)				
	Observed		Standard		
HPLC Purity:	99% at 254 nm		≥ 90% at 254 nm		
Absorption (Buffer pH 10)					
Wavelength (λ _{max}):	373 nm		372 ± 3 nm		
Extinction:	17,400 cm ⁻¹ M ⁻¹		16,500 ± 2,000 cm ⁻¹ M ⁻¹		
TLC Purity:	> 99%		≥ 90%		
TLC System:	1:1 Hexane/ Ethyl acetate				
Storage:	Refrigerate: Store at 4-8°C (Ref), Protect from moisture (D), Protect from light (L)				
Production Lot:	14527	Release:	03-Nov-2014	Expiration:	03-Nov-2024
Structure:					

This data reflects the accuracy of the analysis of the subject lot
 For research use only. Not for human use. No animal source.
 This material is chemically synthesized and manufactured in the USA

Approved by: Dianna M Wineinger
 Analytical Services Manager

Date: 03-Nov-2014

Signature: 