


Certificate of Analysis (C of A)								
Name:		EverFluor FL Cholesterol						
Catalog Number:		7356	Lot:		50241			
CAS Number:		N/A	Formula:		C <sub>36</sub> H <sub>50</sub> BF <sub>2</sub> N <sub>3</sub> O <sub>2</sub>			
Molecular Weight:		605.62 (g/mol)						
		Observed		Standard				
HPLC Purity:		99% at 503 nm		≥ 95% at 503 nm				
Absorption (Methanol)								
Wavelength (λ <sub>max</sub> ):		504 nm		503 ± 3 nm				
Extinction:		Set 88,000 cm <sup>-1</sup> M <sup>-1</sup>		Set 88,000 ± N/A cm <sup>-1</sup> M <sup>-1</sup>				
TLC Purity:		> 99%		≥ 95%				
TLC System:		90:10 Chloroform / Methanol						
Comments:		Set ABS Extinction at 88,000.						
Storage:		Freeze: Store at -5 to -30°C (F), Protect from moisture (D), Protect from light (L)						
Production Lot:		14652	Release:		04-Dec-2014	Expiration:		04-Dec-2024
Structure:								

This data reflects the accuracy of the analysis of the subject lot  
 For research use only. Not for human use. No animal source.  
 This material is chemically synthesized and manufactured in the USA

**Approved by:** Dianna M Wineinger  
 Analytical Services Manager

**Date:** 04-Dec-2014

**Signature:**

*Dianna M Wineinger*