


Certificate of Analysis (C of A)					
Name:	mBCI Monochlorobimane				
Catalog Number:	6063		Lot:	51102	
CAS Number:	76421-73-3		Formula:	C ₁₀ H ₁₁ ClN ₂ O ₂	
Molecular Weight:	226.66 (g/mol)				
	Observed	Standard			
HPLC Purity:	98% at 254 nm		≥ 95% at 254 nm		
Absorption (Methanol)					
Wavelength (λ_{max}):	379 nm		380 ± 3 nm		
Extinction:	6,100 cm ⁻¹ M ⁻¹		6,200 ± 600 cm ⁻¹ M ⁻¹		
TLC Purity:	> 99%		≥ 95%		
TLC System:	90:10 Chloroform/Methanol				
Storage:	Freeze: Store at -5 to -30°C (F), Protect from moisture (D), Protect from light (L)				
Production Lot:	15695	Release:	16-Feb-2017	Expiration:	16-Feb-2027
Structure:					

This data reflects the accuracy of the analysis of the subject lot
 For research use only. Not for human use. No animal source.
 This material is chemically synthesized and manufactured in the USA

Approved by: Dianna M Wineinger
 Analytical Services Manager

Date: 16-Feb-2017

Signature: 