


Certificate of Analysis (C of A)					
Name:		CellHunt Calcein Green, AM *Sp-Pak*			
Catalog Number:		7912	Lot:		51210
CAS Number:		148504-34-1	Formula:		C ₄₆ H ₄₆ N ₂ O ₂₃
Molecular Weight:		994.87 (g/mol)			
		Observed	Standard		
HPLC Purity:		96% at 254 nm	≥ 80% at 254 nm		
Absorption (Methanol)					
Wavelength (λ _{max}):		283 nm	283 ± 3 nm		
Extinction:		5,100 cm ⁻¹ M ⁻¹	5,500 ± 500 cm ⁻¹ M ⁻¹		
TLC Purity:		> 98%	≥ 80%		
TLC System:		60:40 Ethyl acetate / Hexane			
Storage:		Freeze: Store at -5 to -30°C (F), Protect from moisture (D), Protect from light (L)			
Production Lot:		14818	Release:	10-Mar-2015	Expiration: 02-Dec-2024
Structure:					

This data reflects the accuracy of the analysis of the subject lot
 For research use only. Not for human use. No animal source.
 This material is chemically synthesized and manufactured in the USA

Approved by: Dianna M Wineinger
 Analytical Services Manager

Date: 10-Mar-2015

Signature: 