


Certificate of Analysis (C of A)			
Name:	Calcein, AM		
Catalog Number:	7077	Lot:	52239
CAS Number:	148504-34-1	Formula:	C ₄₆ H ₄₆ N ₂ O ₂₃
Molecular Weight:	994.87 (g/mol)	Appearance:	White Solid
	Observed	Standard	
HPLC Purity:	94% at 254 nm	≥ 80% at 254 nm	
Absorption (Methanol)			
Wavelength (λ_{max}):	283 nm	283 ± 3 nm	
Extinction:	5,100 cm ⁻¹ M ⁻¹	5,500 ± 500 cm ⁻¹ M ⁻¹	
TLC Purity:	> 98%	≥ 80%	
TLC System:	60:40 Ethyl acetate / Hexane		
Storage:	Freeze: Store at -5 to -30°C (F), Protect from moisture (D), Protect from light (L)		
Production Lot:	14643	Release:	06-Apr-2020
		Expiration:	02-Dec-2024
Structure:			

This data reflects the accuracy of the analysis of the subject lot
 For research use only. Not for human use. No animal source.
 This material is chemically synthesized and manufactured in the USA

Approved by: Dianna M Wineinger
 Analytical Services Manager

Date: 06-Apr-2020

Signature: 