


Certificate of Analysis (C of A)					
Name:		CellHunt Calcein Violet, AM for 405 nm excitation			
Catalog Number:	7612	Lot:	50834		
Molecular Weight:	616.48 (g/mol)				
		Observed		Standard	
HPLC Purity:	95% at 301 nm		≥ 90% at 301 nm		
Absorption (Acetonitrile)					
Wavelength (λ_{\max}):	297 nm		298 ± 3 nm		
Extinction:	Set 13,000 cm ⁻¹ M ⁻¹		Set 13,000 ± 1,000 cm ⁻¹ M ⁻¹		
Storage:	Freeze: Store at -5 to -30°C (F), Protect from moisture (D), Protect from light (L)				
Production Lot:	15351	Release:	30-Mar-2020	Expiration:	30-Mar-2025

This data reflects the accuracy of the analysis of the subject lot
 For research use only. Not for human use. No animal source.
 This material is chemically synthesized and manufactured in the USA

Approved by: Dianna M Wineinger
 Analytical Services Manager

Date: 30-Mar-2020

Signature: 

Date: 26-Sep-2022