



Certificate of Analysis (C of A)					
Name:	AMS 4-Acetamido-4'-maleimidylstilbene-2,2'-disulfonic acid, disodium salt				
Catalog Number:	6508		Lot:	50359	
CAS Number:	175171-43-4		Formula:	C ₂₀ H ₁₄ N ₂ Na ₂ O ₉ S ₂	
Molecular Weight:	536.44 (g/mol)				
Observed vs Standard					
HPLC Purity:	94% at 254 nm		≥ 65% at 254 nm		
Absorption (Buffer pH 8)					
Wavelength (λ_{max}):	322 nm		323 ± 3 nm		
Extinction:	33,800 cm ⁻¹ M ⁻¹		31,500 ± 4,500 cm ⁻¹ M ⁻¹		
TLC Purity:	> 95%		≥ 65%		
TLC System:	45:45:10 Chloroform/ Methanol/ Acetic acid				
Storage:	Freeze: Store at -5 to -30°C (F), Protect from moisture (D), Protect from light (L)				
Production Lot:	14519	Release:	28-Oct-2014	Expiration:	28-Oct-2024
Structure:					

This data reflects the accuracy of the analysis of the subject lot
 For research use only. Not for human use. No animal source.
 This material is chemically synthesized and manufactured in the USA

Approved by: Dianna M Wineinger
 Analytical Services Manager

Date: 28-Oct-2014

Signature: 

Date: 09-Dec-2022