


Certificate of Analysis (C of A)					
Name:	D-Biotin, SE D-Biotin, succinimidyl ester				
Catalog Number:	6739		Lot:	50780	
CAS Number:	35013-72-0		Formula:	C ₁₄ H ₁₉ N ₃ O ₅ S	
Molecular Weight:	341.38 (g/mol)				
	Observed		Standard		
TLC Purity:	> 99%		≥ 95%		
TLC System:	90:10 Chloroform / Methanol with Biotin Dip				
Storage:	Freeze: Store at -5 to -30°C (F), Protect from moisture (D), Protect from light (L)				
Production Lot:	15325	Release:	17-Mar-2016	Expiration:	17-Mar-2021
Structure:					

This data reflects the accuracy of the analysis of the subject lot
 For research use only. Not for human use. No animal source.
 This material is chemically synthesized and manufactured in the USA

Approved by: Dianna M Wineinger
 Analytical Services Manager

Date: 17-Mar-2016

Signature: 