


Certificate of Analysis (C of A)					
Name:	Biocytin ε-Biotinoyl-L-lysine				
Catalog Number:	6718		Lot:	51793	
CAS Number:	576-19-2		Formula:	C ₁₆ H ₂₈ N ₄ O ₄ S	
Molecular Weight:	372.48 (g/mol)				
	Observed			Standard	
TLC Purity:	> 99%			≥ 95%	
TLC System:	80:20 Acetonitrile / Water with Biotin Dip and Ninhydrin Dip				
Storage:	Refrigerate: Store at 4-8°C (Ref), Protect from moisture (D), Protect from light (L)				
Production Lot:	15073	Release:	11-Aug-2015	Expiration:	11-Aug-2025
Structure:					

This data reflects the accuracy of the analysis of the subject lot
 For research use only. Not for human use. No animal source.
 This material is chemically synthesized and manufactured in the USA

Approved by: Dianna M Wineinger
 Analytical Services Manager

Date: 11-Aug-2015

Signature: 