


Certificate of Analysis (C of A)					
Name:	CLZN hcp Coelenterazine hcp				
Catalog Number:	6760		Lot:	50288	
CAS Number:	123437-32-1		Formula:	C ₂₅ H ₂₅ N ₃ O ₂	
Molecular Weight:	399.48 (g/mol)				
	Observed		Standard		
HPLC Purity:	97% at 254 nm		≥ 90% at 254 nm		
Absorption (Methanol)					
Wavelength (λ_{max}):	429 nm		430 ± 3 nm		
Extinction:	8,200 cm ⁻¹ M ⁻¹		8,500 ± 800 cm ⁻¹ M ⁻¹		
TLC Purity:	> 98%		≥ 90%		
TLC System:	90:10 Chloroform / Methanol				
Storage:	Freeze: Store at -5 to -30°C (F), Protect from moisture (D), Protect from light (L)				
Production Lot:	14728	Release:	22-Jan-2015	Expiration:	22-Jan-2020
Structure:					

This data reflects the accuracy of the analysis of the subject lot
 For research use only. Not for human use. No animal source.
 This material is chemically synthesized and manufactured in the USA

Approved by: Dianna M Wineinger
 Analytical Services Manager

Date: 22-Jan-2015

Signature: 