


| Certificate of Analysis (C of A) | | | | | |
|-------------------------------------|---|----------|-------------|---|-------------|
| Name: | 3β-Amino-3-deoxydigoxigenin hemisuccinamide, succinimidyl ester | | | | |
| Catalog Number: | 6970 | | Lot: | 52138 | |
| CAS Number: | 216299-46-6 | | Formula: | C ₃₁ H ₄₂ N ₂ O ₉ | |
| Molecular Weight: | 586.67 (g/mol) | | | | |
| | | | | | |
| | Observed | | | Standard | |
| TLC Purity: | > 95% | | | ≥ 80% | |
| TLC System: | 90:10 Chloroform / Methanol + heat | | | | |
| Comments: | TLC: Visualized with heat | | | | |
| Storage: | Freeze: Store at -5 to -30°C (F), Protect from moisture (D), Protect from light (L) | | | | |
| Production Lot: | 15311 | Release: | 24-Feb-2021 | Expiration: | 24-Feb-2026 |
| Structure: |  | | | | |

This data reflects the accuracy of the analysis of the subject lot
 For research use only. Not for human use. No animal source.
 This material is chemically synthesized and manufactured in the USA

Approved by: Dianna M Wineinger
 Analytical Services Manager

Date: 24-Feb-2021

Signature: 