


Certificate of Analysis (C of A)					
Name:	PPT 97 β-D-Galactoside PPT 97 β-D-Galactopyranoside				
Catalog Number:	7350	Lot:	51011		
CAS Number:	N/A	Formula:	C ₂₀ H ₁₈ Cl ₂ N ₂ O ₇		
Molecular Weight:	469.28 (g/mol)				
	Observed	Standard			
HPLC Purity:	98% at 254 nm	≥ 90% at 254 nm			
Absorption (0.2M KOH/Methanol)					
Wavelength (λ_{max}):	298 nm	298 ± 3 nm			
Extinction:	N/A cm ⁻¹ M ⁻¹	N/A ± cm ⁻¹ M ⁻¹			
Storage:	Freeze: Store at -5 to -30°C (F), Protect from moisture (D), Protect from light (L)				
Production Lot:	15469	Release:	11-Nov-2016	Expiration:	11-Nov-2026
Structure:					

This data reflects the accuracy of the analysis of the subject lot
 For research use only. Not for human use. No animal source.
 This material is chemically synthesized and manufactured in the USA

Approved by: Dianna M Wineinger
 Analytical Services Manager

Date: 11-Nov-2016

Signature:

Dianna M. Wineinger