


Certificate of Analysis (C of A)					
Name:	AMS 4-Acetamido-4'-maleimidylstilbene-2,2'-disulfonic acid, disodium salt				
Catalog Number:	6508		Lot:	52347	
CAS Number:	175171-43-4		Formula:	C ₂₀ H ₁₄ N ₂ Na ₂ O ₉ S ₂	
Molecular Weight:	536.44 (g/mol)				
	Observed		Standard		
HPLC Purity:	89% at 254 nm		≥ 65% at 254 nm		
Absorption (Buffer pH8)					
Wavelength (λ _{max}):	322 nm		323 ± 3 nm		
Extinction:	29,200 cm ⁻¹ M ⁻¹		31,500 ± 4,500 cm ⁻¹ M ⁻¹		
TLC Purity:	> 98%		≥ 65%		
TLC System:	45:45:10 Chloroform / Methanol / Acetic acid				
Storage:	Freeze: Store at -5 to -30°C (F), Protect from moisture (D), Protect from light (L)				
Production Lot:	15919	Release:	01-Nov-2017	Expiration:	01-Nov-2027
Structure:					

This data reflects the accuracy of the analysis of the subject lot
 For research use only. Not for human use. No animal source.
 This material is chemically synthesized and manufactured in the USA

Approved by: Dianna M Wineinger
 Analytical Services Manager

Date: 01-Nov-2017

Signature: Dianna M Wineinger