


Certificate of Analysis (C of A)					
Name:	Biocytin hydrazide Link-EZ Hydrazide-Biocytin				
Catalog Number:	6769	Lot:	50176		
CAS Number:	102743-85-1	Formula:	C ₁₆ H ₃₀ N ₆ O ₃ S		
Molecular Weight:	386.51 (g/mol)				
	Observed	Standard			
TLC Purity:	> 97%	≥ 90%			
TLC System:	60:40 Acetonitrile / Water with Biotin Dip				
Storage:	Refrigerate: Store at 4-8°C (Ref), Protect from moisture (D), Protect from light (L)				
Production Lot:	11929	Release:	09-Jul-2020	Expiration:	24-Sep-2024
Structure:					

This data reflects the accuracy of the analysis of the subject lot
 For research use only. Not for human use. No animal source.
 This material is chemically synthesized and manufactured in the USA

Approved by: Dianna M Wineinger
 Analytical Services Manager

Date: 09-Jul-2020

Signature: 