


Certificate of Analysis (C of A)

Name:		3β-Amino-3-deoxydigoxigenin hemisuccinamide, succinimidyl ester				
Catalog Number:		6970		Lot: 52138		
CAS Number:		216299-46-6		Formula: C ₃₁ H ₄₂ N ₂ O ₉		
Molecular Weight:		586.67 (g/mol)				
		Observed		Standard		
TLC Purity:		> 95%		≥ 80%		
TLC System:		90:10 Chloroform / Methanol + heat				
Comments:		TLC: Visualized with heat				
Storage:		Freeze: Store at -5 to -30°C (F), Protect from moisture (D), Protect from light (L)				
Production Lot:		15311	Release:	24-Feb-2021	Expiration:	24-Feb-2026
Structure:						

This data reflects the accuracy of the analysis of the subject lot
For research use only. Not for human use. No animal source.
This material is chemically synthesized and manufactured in the USA

Approved by: Dianna M Wineinger
Analytical Services Manager

Date: 24-Feb-2021

Signature:

Dianna M Wineinger