



Certificate of Analysis (C of A)					
Name:	N-(1-Pyrene)iodoacetamide Acetamide, 2-iodo-N-(1-pyrenyl)				
Catalog Number:	6947		Lot:	50114	
CAS Number:	76936-87-3		Formula:	C <sub>18</sub> H <sub>12</sub> INO	
Molecular Weight:	385.20 (g/mol)				
	<b>Observed</b>		<b>Standard</b>		
HPLC Purity:	93% at 254 nm		≥ 90% at 254 nm		
<b>Absorption (Methanol)</b>					
Wavelength (λ <sub>max</sub> ):	340 nm		339 ± 3 nm		
Extinction:	29,400 cm <sup>-1</sup> M <sup>-1</sup>		27,500 ± 2,000 cm <sup>-1</sup> M <sup>-1</sup>		
<b>Fluorescence (Methanol)</b>					
Emission Max:	384 nm		385 ± 4 nm		
Storage:	Freeze: Store at -5 to -30°C (F), Protect from moisture (D), Protect from light (L)				
Production Lot:	14404	Release:	04-Aug-2014	Expiration:	04-Aug-2024
Structure:					

This data reflects the accuracy of the analysis of the subject lot  
 For research use only. Not for human use. No animal source.  
 This material is chemically synthesized and manufactured in the USA

**Approved by:** Nabi  
 Analytical Services Manager  
**Signature:** 

**Date:** 04-Aug-2014