


Certificate of Analysis (C of A)					
Name:	5-IAF 5-Iodoacetamidofluorescein				
Catalog Number:	6881		Lot:	50070	
CAS Number:	63368-54-7		Formula:	C ₂₂ H ₁₄ INO ₆	
Molecular Weight:	515.25 (g/mol)				
	Observed	Standard			
HPLC Purity:	% at nm	≥ 90% at 254 nm			
Absorption ()					
Wavelength (λ _{max}):	0 nm	492 ± 3 nm			
Extinction:	N/A cm ⁻¹ M ⁻¹	78,000 ± 5,000 cm ⁻¹ M ⁻¹			
Fluorescence ()					
Emission Max:		515 ± 4 nm			
Storage:	Freeze: Store at -5 to -30°C (F), Protect from moisture (D), Protect from light (L)				
Production Lot:	ST11420	Release:	07-Feb-2013	Expiration:	07-Feb-2023
Structure:					

This data reflects the accuracy of the analysis of the subject lot
 For research use only. Not for human use. No animal source.
 This material is chemically synthesized and manufactured in the USA

Approved by:

Date: 07-Feb-2013

Analytical Services Manager

Signature: 