


Certificate of Analysis

(C of A)

Name:		NIETf 2-(2,3-Naphthalimino)ethyl trifluoromethanesulfonate						
Catalog Number:		6369	Lot:		52190			
CAS Number:		128651-50-3	Formula:		C ₁₅ H ₁₀ F ₃ NO ₅ S			
Molecular Weight:		373.30 (g/mol)	Appearance:		Off white Solid			
		Observed		Standard				
Absorption (5:1 Methanol / Water)								
Wavelength (λ _{max}):		260 nm		260 ± 3 nm				
Extinction:		60,100 cm ⁻¹ M ⁻¹		63,000 ± 5,000 cm ⁻¹ M ⁻¹				
TLC Purity:		> 98%		≥ 90%				
TLC System:		90:10 Chloroform / Hexane						
Storage:		Freeze: Store at -5 to -30°C (F), Protect from moisture (D), Protect from light (L)						
Production Lot:		16505	Release:		21-Feb-2022	Expiration:		21-Feb-2032
Structure:								

This data reflects the accuracy of the analysis of the subject lot
 For research use only. Not for human use. No animal source.
 This material is chemically synthesized and manufactured in the USA

Approved by: Dianna M. Wineinger
 Analytical Services Manager

Date: 21-Feb-2022

Signature: 