


Certificate of Analysis (C of A)					
Name:		DDAOG DDAO galactoside			
Catalog Number:		6899	Lot:		50544
CAS Number:		503178-95-8	Formula:		C ₂₁ H ₂₁ Cl ₂ NO ₇
Molecular Weight:		470.31 (g/mol)			
		Observed		Standard	
HPLC Purity:		97% at 254 nm		≥ 90% at 254 nm	
Absorption (Buffer pH7)					
Wavelength (λ _{max}):		465 nm		465 ± 3 nm	
Extinction:		27,900 cm ⁻¹ M ⁻¹		29,000 ± 2,500 cm ⁻¹ M ⁻¹	
TLC Purity:		> 99%		≥ 90%	
TLC System:		80:20 Chloroform / Methanol			
Storage:		Freeze: Store at -5 to -30°C (F), Protect from moisture (D), Protect from light (L)			
Production Lot:		11513	Release:	07-Aug-2014	Expiration: 17-Apr-2023
Structure:					

This data reflects the accuracy of the analysis of the subject lot
 For research use only. Not for human use. No animal source.
 This material is chemically synthesized and manufactured in the USA

Approved by: Dianna M Wineinger
 Analytical Services Manager

Date: 07-Aug-2014

Signature: 