


Certificate of Analysis

(C of A)

Name:		MGITC Malachite Green Isothiocyanate	
Catalog Number:	6591	Lot:	51843
CAS Number:	147492-82-8	Formula:	C ₂₄ H ₂₄ ClN ₃ O ₄ S
Molecular Weight:	485.98 (g/mol)		
	Observed	Standard	
HPLC Purity:	98% at 254 nm	≥ 90% at 254 nm	
Absorption (Acetonitrile)			
Wavelength (λ _{max}):	629 nm	629 ± 3 nm	
Extinction:	89,200 cm ⁻¹ M ⁻¹	84,000 ± 7,000 cm ⁻¹ M ⁻¹	
TLC Purity:	> 99%	≥ 90%	
TLC System:	80:20 Chloroform / Methanol		
Storage:	Freeze: Store at -5 to -30°C (F), Protect from moisture (D), Protect from light (L)		
Production Lot:	16334	Release:	13-Jul-2020
		Expiration:	13-Jul-2030
Structure:			

This data reflects the accuracy of the analysis of the subject lot
 For research use only. Not for human use. No animal source.
 This material is chemically synthesized and manufactured in the USA

Approved by: Dianna M Wineinger
 Analytical Services Manager

Date: 13-Jul-2020

Signature: 