


<b>Certificate of Analysis</b> (C of A)					
Name:	<b>N-(2-Aminoethyl)maleimide, trifluoroacetic acid salt</b> <b>C2-Maleimide, TFA</b>				
Catalog Number:	7018		Lot:	50721	
CAS Number:	146474-00-2		Formula:	C <sub>8</sub> H <sub>9</sub> F <sub>3</sub> N <sub>2</sub> O <sub>4</sub>	
Molecular Weight:	254.16 (g/mol)				
	<b>Observed</b>		<b>Standard</b>		
TLC Purity:	> 98%		≥ 90%		
TLC System:	60:30:7:3 Chloroform / Methanol / Water / Acetic Acid with Ninhydrin Dip				
Storage:	Freeze: Store at -5 to -30°C (F), Protect from moisture (D), Protect from light (L)				
Production Lot:	15265	Release:	15-Jan-2016	Expiration:	15-Jan-2026
Structure:					

This data reflects the accuracy of the analysis of the subject lot  
 For research use only. Not for human use. No animal source.  
 This material is chemically synthesized and manufactured in the USA

**Approved by:** Dianna M Wineinger  
 Analytical Services Manager

**Date:** 15-Jan-2016

**Signature:** 