


Certificate of Analysis (C of A)			
Name:	PPT 97 β -D-Galactoside PPT 97 β -D-Galactopyranoside		
Catalog Number:	7350	Lot:	51011
CAS Number:	N/A	Formula:	$C_{20}H_{18}Cl_2N_2O_7$
Molecular Weight:	469.28 (g/mol)		
	Observed	Standard	
HPLC Purity:	98% at 254 nm	$\geq 90\%$ at 254 nm	
Absorption (0.2M KOH/Methanol)			
Wavelength (λ_{max}):	298 nm	298 $\pm$ 3 nm	
Extinction:	N/A $cm^{-1}M^{-1}$	N/A $\pm$ $cm^{-1}M^{-1}$	
Storage:	Freeze: Store at -5 to -30°C (F), Protect from moisture (D), Protect from light (L)		
Production Lot:	15469	Release:	11-Nov-2016
		Expiration:	11-Nov-2026
Structure:			

This data reflects the accuracy of the analysis of the subject lot
 For research use only. Not for human use. No animal source.
 This material is chemically synthesized and manufactured in the USA

Approved by: Dianna M Wineinger
 Analytical Services Manager

Date: 11-Nov-2016

Signature:

Dianna M. Wineinger