



Certificate of Analysis (C of A)					
Name:	PBFI, AM *Cell permeant* Potassium-binding benzofuran isophthalate acetoxymethyl ester				
Catalog Number:	6074		Lot:	52360	
CAS Number:	124549-23-1		Formula:	$C_{58}H_{62}N_2O_{24}$	
Molecular Weight:	1171.11 (g/mol)				
	Observed		Standard		
HPLC Purity:	96% at 254 nm		≥ 95% at 254 nm		
Absorption (Methanol)					
Wavelength (λ_{max}):	367 nm		367 ± 3 nm		
Extinction:	Set 36,000 $cm^{-1}M^{-1}$		Set 36,000 ± 3,500 $cm^{-1}M^{-1}$		
TLC Purity:	> 97%		≥ 95%		
TLC System:	90:10 Chlororm / Methanol				
Storage:	Freeze: Store at -5 to -30°C (F), Protect from moisture (D), Protect from light (L)				
Production Lot:	16594	Release:	13-Apr-2023	Expiration:	13-Apr-2033
Structure:					

This data reflects the accuracy of the analysis of the subject lot
 For research use only. Not for human use. No animal source.
 This material is chemically synthesized and manufactured in the USA

Approved by: Nabi
 Analytical Services Manager
Signature: 

Date: 13-Apr-2023