


Certificate of Analysis (C of A)					
Name:	BF 405, SE Bella Fluor 405 NHS Ester				
Catalog Number:	7654	Lot:	52040		
CAS Number:	N/A	Formula:	C ₄₆ H ₆₉ N ₅ O ₁₅ S ₃		
Molecular Weight:	1028.28 (g/mol)				
	Observed	Standard			
HPLC Purity:	93% at 254 nm	≥ 80% at 254 nm			
Absorption (Methanol)					
Wavelength (λ_{max}):	396 nm	397 ± 3 nm			
Extinction:	37,900 cm ⁻¹ M ⁻¹	38,000 ± 2,500 cm ⁻¹ M ⁻¹			
Fluorescence (Methanol)					
Emission Max:	406 nm	406 ± 4 nm			
Storage:	Freeze: Store at -5 to -30°C (F), Protect from moisture (D), Protect from light (L)				
Production Lot:	15123	Release:	03-Mar-2021	Expiration:	03-Mar-2026
Structure:	 <p style="text-align: center;">3 (CH₃CH₂)₃NH</p>				

This data reflects the accuracy of the analysis of the subject lot
 For research use only. Not for human use. No animal source.
 This material is chemically synthesized and manufactured in the USA

Approved by: Dianna M Wineinger
 Analytical Services Manager

Date: 03-Mar-2021

Signature: 